

Anti-Phospho-ErbB3/Her3-Y1328 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21860

Item Number: CM0001355



Introduction

Phospho-specific rabbit polyclonal antibody targeting ErbB3/HER3 phosphorylated at Y1328, validated for WB, IF/ICC, and ELISA in human, mouse, and rat. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ErbB3/HER3-Y1328 Rabbit pAb
Remarks/Alias	HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; Phospho-ErbB3/HER3-Y1328
Species	Human
Gene ID (Human)	2065
Gene ID	2065
Immunogen	A synthetic phosphorylated peptide around Y1328 of human ERBB3 (NP_001973.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21860
Protein Weight	148kDa
Shipping	Ice bag