


Anti-Phospho-ErbB4/Her4-Y1284 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15303

Item Number: CM0001819



Introduction

Phospho-specific rabbit monoclonal antibody against ErbB4/HER4 phosphorylated at tyrosine 1284. Validated for WB, IF/ICC, and ELISA in human samples; ideal for studying ErbB4-mediated signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ErbB4/HER4-Y1284 Rabbit mAb
Remarks/Alias	HER4; ALS19; p180erbB4
Species	Human
Gene ID	2066
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y1284 of human ErbB4/HER4 (NP_005226.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15303
Protein Weight	147 kDa
Shipping	Ice bag