


Anti-Phospho-Erk1 (T202) + Erk2 (T185) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P28482 / P27361

Item Number: CM0001133



Introduction

This phospho-specific rabbit monoclonal antibody detects ERK1 phosphorylated at T202 and ERK2 at T185. Validated for WB, IHC-P, and ELISA in human, mouse, and rat. Ideal for MAPK/ERK signaling pathway studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ERK1-T202 + ERK2-T185 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; H3/A; H3FA; ERK1 / ERK2; Phospho-ERK1-T202 + ERK2-T185
Species	Human
Gene ID (Human)	5594/5595
Gene ID	5594/5595
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T202 of human ERK1 (P27361).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	42,44kDa
Shipping	Ice bag