


Anti-Phospho-Erk1-T202/Y204 + Erk2-T185/Y187 Rabbit Pab For Wb, If, Elisa - P28482 / P27361

Item Number: CM0001108



Introduction

Phospho-ERK1 (T202/Y204) + ERK2 (T185/Y187) rabbit polyclonal antibody for WB, IF/ICC, ELISA. Human, mouse, rat reactivity. Ideal for MAPK/ERK signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-ERK1-T202/Y204 + ERK2-T185/Y187 Rabbit pAb
Remarks/Alias	MAPK1/MAPK3; Phospho-ERK1-T202/Y204 + ERK2-T185/Y187
Species	Human
Gene ID (Human)	5594/5595
Gene ID	5594/5595
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28482/P27361
Protein Weight	36kDa/41kDa/38kDa/40kDa/43kDa
Shipping	Ice bag