


Anti-Phospho-Esr α -S118 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P03372

Item Number: CM0000770



Introduction

Phospho-ESR α -S118 Rabbit Polyclonal Antibody detects phosphorylated Estrogen Receptor alpha (Ser118) in human samples. Applications: WB, IF/ICC, ELISA. Molecular weight: 66kDa.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ESR α -S118 Rabbit pAb
Remarks/Alias	ER; ESR; Era; ESRA; ESTRR; NR3A1; Phospho-ESR α -S118
Species	Human
GeneID (human)	2099
GeneID	2099
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around S118 of human Estrogen Receptor alpha (ESR1) (NP_000116.2).)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P03372
Protein Weight	66kDa
Shipping	Ice bag