


Anti-Phospho-Fak-Y576/Y577 Rabbit Monoclonal Antibody For Wb, Elisa - Q05397

Item Number: CM0001390



Introduction

Phospho-specific rabbit monoclonal antibody detecting FAK phosphorylated at Y576 and Y577. Validated for Western blotting and ELISA. Cross-reacts with mouse. Essential for research on focal adhesion signaling, cell migration, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FAK-Y576/Y577 Rabbit mAb
Remarks/Alias	FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK
Species	Human
Gene ID (Human)	5747
Gene ID	5747
Immunogen	A synthetic phosphorylated peptide around Y576 & Y577 of human FAK (NP_722560.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q05397
Protein Weight	119kDa
Shipping	Ice bag