


## Anti-Phospho-Fak-Y925 Rabbit Polyclonal Antibody For Wb, Elisa - Q05397

Item Number: CM0001512



### Introduction

Phospho-specific rabbit pAb targeting FAK (Y925). Validated for human samples in WB and ELISA. Detects focal adhesion kinase 1 phosphorylated at tyrosine 925, essential for cell migration, adhesion, and integrin signaling research. SWISS: Q05397.

[Learn More](#)

| Parameter        | Value                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------|
| Product Name     | Phospho-FAK-Y925 Rabbit pAb                                                                   |
| Remarks/Alias    | FAK; FADK; FAK1; FRNK; FADK 1; PPP1R71; p125FAK; pp125FAK; Phospho-FAK-Y925                   |
| Species          | Human                                                                                         |
| Gene ID (Human)  | 5747                                                                                          |
| Gene ID          | 5747                                                                                          |
| Immunogen        | Synthetic peptide, A synthetic phosphorylated peptide around Y925 of human FAK (NP_722560.1). |
| Source           | Rabbit                                                                                        |
| Category         | Phosphorylated Antibodies                                                                     |
| Application      | WB, ELISA                                                                                     |
| Cross Reactivity | Human                                                                                         |
| SWISS            | Q05397                                                                                        |
| Protein Weight   | 119kDa                                                                                        |
| Shipping         | Ice bag                                                                                       |