

Anti-Phospho-Fgfr1 (Y653) Rabbit Polyclonal Antibody For Wb, Elisa - P11362

Item Number: CM0000613



Introduction

Phospho-FGFR1 (Y653) Rabbit Polyclonal Antibody validated for WB, ELISA in Human, Mouse, Rat. Detects phosphorylated FGFR1 at tyrosine 653, a receptor tyrosine kinase essential for embryonic development and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FGFR1-Y653 Rabbit pAb
Remarks/Alias	CEK; FLG; HH2; OGD; ECCL; FLT2; KAL2; BFGFR; CD331; FGFBR; FLT-2; HBGFR; N-SAM; FGFR-1; HRTFDS; bFGF-R-1; Phospho-FGFR1-Y653
Species	Human
GeneID (Human)	2260
GeneID	2260
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y653 of human FGFR1 (NP_075598.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11362
Protein Weight	92kDa
Shipping	Ice bag