

Anti-Phospho-Fgfr1 (Y653/Y654) Rabbit Polyclonal Antibody For Wb, Elisa - P11362

Item Number: CM0001065



Introduction

Rabbit polyclonal antibody specific for phosphorylated FGFR1 (Tyr653/Tyr654). Verified for Western blot and ELISA in human and mouse samples. Ideal for studying FGFR signaling in development and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FGFR1-Y653/Y654 Rabbit pAb
Remarks/Alias	CEK; FLG; HH2; OGD; ECCL; FLT2; KAL2; BFGFR; CD331; FGFBR; FLT-2; HBGFR; N-SAM; FGFR-1; HRTFDS; bFGF-R-1; Phospho-FGFR1-Y653/Y654
Species	Human
GeneID (Human)	2260
GeneID	2260
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y653 & Y654 of human Phospho-FGFR1-Y653/Y654 (NP_075598.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P11362
Protein Weight	92kDa
Shipping	Ice bag