


Anti-Phospho-Flt3-Y589/591 Rabbit Polyclonal Antibody For Wb, Elisa - P36888

Item Number: CM0001666



Introduction

Phospho-FLT3-Y589/591 Rabbit Polyclonal Antibody specifically detects phosphorylation at Y589/591 of human FLT3 receptor tyrosine kinase. Suitable for Western blot and ELISA, it supports research into hematopoietic signaling and acute myeloid leukemia.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FLT3-Y589/591 Rabbit pAb
Remarks/Alias	FLK2; STK1; CD135; FLK-2; Phospho-FLT3-Y589/591
Species	Human
Gene ID (Human)	2322
Gene ID	2322
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y589 & Y591 of human FLT3 (NP_004110.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P36888
Protein Weight	113kDa
Shipping	Ice bag