

Anti-Phospho-Foxo1-T24/Foxo3A-T32/Foxo4-T28 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q12778 / O43524 / P98177

Item Number: CM0001794



Introduction

Phospho-specific rabbit polyclonal antibody recognizing FOXO1 (T24), FOXO3A (T32), and FOXO4 (T28). Ideal for studying insulin signaling, apoptosis, and metabolic regulation in human, mouse, and rat samples via WB, IHC-P, and ELISA.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FOXO1-T24/FOXO3A-T32/FOXO4-T28 Rabbit pAb
Remarks/Alias	FoxO1/FoxO3a/FoxO4; Phospho-FOXO1-T24/FOXO3A-T32/FOXO4-T28
Species	Human
Gene ID (Human)	2308/2309/4303
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T24 of human FOXO1/FOXO3A/FOXO4 (FOXO1 NP_002006.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12778/O43524/P98177
Shipping	Ice bag