

Anti-Phospho-Foxo3A (Ser294) Polyclonal Antibody For Wb, Elisa - O43524

Item Number: CM0001537



Introduction

Phospho-FOXO3A (Ser294) Rabbit pAb detects human, mouse, rat FOXO3A phosphorylated at S294. Suitable for WB and ELISA. Ideal for studying apoptosis, autophagy, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-FOXO3A-S294 Rabbit pAb
Remarks/Alias	FOXO2; AF6q21; FKHL1; FOXO3A; FKHL1P2; Phospho-FOXO3A-S294
Species	Human
Gene ID (Human)	2309
Gene ID	2309
Immunogen	A synthetic phosphorylated peptide around S294 of human FOXO3A (NP_001446.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O43524
Protein Weight	71kDa
Shipping	Ice bag