

Anti-Phospho-Fyn-S21 Rabbit Polyclonal Antibody For Wb, Elisa - P06241

Item Number: CM0001023



Introduction

Rabbit polyclonal antibody targeting Phospho-Fyn-S21 (human, 61kDa) for WB and ELISA. Ideal for researching Fyn kinase roles in cell adhesion, immune response, and neural signaling.

[Learn More](#)

Specification	Details
Product Name	Phospho-Fyn-S21 Rabbit pAb
Remarks/Alias	SLK; SYN; p59-FYN; Phospho-Fyn-S21
Species	Human
Gene ID (Human)	2534
Gene ID	2534
Immunogen	A synthetic phosphorylated peptide around S21 of human Phospho-Fyn-S21 (NP_002028.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06241
Protein Weight	61kDa
Shipping	Ice bag