


Anti-Phospho-Gab1 (Y307) Rabbit Polyclonal Antibody For Wb, Elisa - Q13480

Item Number: CM0000762



Introduction

Rabbit polyclonal antibody targeting phosphorylated Y307 of human GAB1, suitable for WB and ELISA. Detects mouse GAB1 (predicted 77 kDa). Ideal for studying EGFR, insulin receptor, and MET/HGF signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GAB1-Y307 Rabbit pAb
Remarks/Alias	DFNB26; Phospho-GAB1-Y307
Species	Human
Gene ID (Human)	2549
Gene ID	2549
Immunogen	A synthetic phosphorylated peptide around Y307 of human GAB1 (NP_002030.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13480
Protein Weight	77kDa
Shipping	Ice bag