


Anti-Phospho-Gcn2 (T899) Rabbit Polyclonal Antibody For Wb, Elisa - Q9P2K8

Item Number: CM0001994



Introduction

Phospho-GCN2-T899 Rabbit Polyclonal Antibody validated for Western blotting and ELISA. Detects human GCN2 phosphorylated at threonine 899. Ideal for integrated stress response and amino acid starvation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GCN2-T899 Rabbit pAb
Remarks/Alias	GCN2; PVOD2; Phospho-GCN2-T899
Species	Human
GeneID (human)	440275
GeneID	440275
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T899 of human Phospho-GCN2-T899.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P2K8
Protein Weight	187kDa
Shipping	Ice bag