


Anti-Phospho-GluR1/Gria1-S849 Rabbit Polyclonal Antibody For Wb, Elisa - P42261

Item Number: CM0001138



Introduction

Rabbit polyclonal antibody targeting phosphorylated Ser849 of human GluR1/GRIA1. Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. An essential tool for studying glutamate receptor signaling in neuroscience.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GluR1/GRIA1-S849 Rabbit pAb
Remarks/Alias	GLUH1; GLUR1; GLURA; GluA1; HBGR1; MRD67; MRT76; Phospho-GluR1/GRIA1-S849
Species	Human
Gene ID (Human)	2890
Gene ID	2890
Immunogen	Synthetic peptide: a phospho specific peptide corresponding to residues surrounding S849 of human GluR1/GRIA1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42261
Protein Weight	102kDa
Shipping	Ice bag