


Anti-Phospho-Gsk3A (S21) Rabbit Polyclonal Antibody - P49840

Item Number: CM0000607



Introduction

High-quality rabbit polyclonal antibody targeting phosphorylated GSK3α at Ser21 (Human) for WB, IP, ELISA. Ideal for studying Wnt signaling, insulin resistance, apoptosis, and Alzheimer's disease. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3α-S21 Rabbit pAb
Remarks/Alias	GSK3 alpha; GSK3A; Phospho-GSK3α-S21
Species	Human
GeneID (Human)	2931
GeneID	2931
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S21 of human GSK3α (NP_063937.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P49840
Protein Weight	51kDa
Shipping	Ice bag