


Anti-Phospho-Gsk3B-S9 Monoclonal Antibody For Wb, Elisa - P49841

Item Number: CM0001502



Introduction

High-quality Phospho-GSK3 β -S9 Mouse Monoclonal Antibody validated for WB and ELISA. Detects human, mouse, and rat GSK3 β phosphorylated at Ser9. Ideal for studying Wnt signaling, glucose metabolism, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GSK3 β -S9 Mouse mAb
Remarks/Alias	GSK3B; gsk-3 β ; Phospho-GSK3 β -S9
Species	Human
Gene ID (Human)	2932
Gene ID	2932
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding S9 of human GSK3 β
Source	Mouse
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag