


Anti-Phospho-GYS1-S641 Rabbit Polyclonal Antibody For Wb, Elisa - P13807

Item Number: CM0001840



Introduction

Rabbit polyclonal antibody targeting GYS1 phosphorylated at Ser641, validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying glycogen synthesis and metabolic regulation. ~84 kDa band detected.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GYS1-S641 Rabbit pAb
Remarks/Alias	GYS; GYS; Phospho-GYS1-S641
Species	Human
GeneID (human)	2997
GeneID	2997
Immunogen	Synthetic phosphorylated peptide around S641 of human Glycogen synthase 1 (GYS1) (NP_002094.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13807
Protein Weight	84kDa
Shipping	Ice bag