


Anti-Phospho-Hdac4-S246 + Hdac5-S259 + Hdac9-S220 Rabbit Monoclonal Antibody - P56524 / Q9Uql6 / Q9Ukv0

Item Number: CM0001682



Introduction

Rabbit monoclonal antibody targeting phosphorylated HDAC4 (S246), HDAC5 (S259), and HDAC9 (S220). Suitable for WB, IHC-P, ELISA in human, mouse, and rat samples. Chromatin regulator and transcription repressor.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC4-S246 + HDAC5-S259 + HDAC9-S220 Rabbit mAb
Remarks/Alias	HDAC4/HDAC5/HDAC9; Phospho-HDAC4-S246 + HDAC5-S259 + HDAC9-S220
Species	Human
GeneID (Human)	9759/10014/9734
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56524/Q9UQL6/Q9UKV0
Protein Weight	57-65 kDa/97-122 kDa
Shipping	Ice bag