

## Anti-Phospho-Her2/ErbB2 (Y1221/1222) Rabbit Polyclonal Antibody For Wb, Elisa - P04626

Item Number: CM0000825



### Introduction

Phospho-specific rabbit polyclonal antibody detecting HER2/ErbB2 phosphorylated at Y1221/1222. Validated for Western blot (WB) and ELISA applications in human samples. Ideal for oncology and signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                             |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-HER2/ErbB2-Y1221/1222 Rabbit pAb                                                                                          |
| Remarks/Alias    | NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2); Phospho-HER2/ErbB2-Y1221/1222 |
| Species          | Human                                                                                                                             |
| GeneID (human)   | 2064                                                                                                                              |
| GeneID           | 2064                                                                                                                              |
| Immunogen        | Synthetic peptide: A synthetic phosphorylated peptide around Y1221 & Y1222 of human HER2/ERBB2 (NP_004439.2)                      |
| Source           | Rabbit                                                                                                                            |
| Category         | Phosphorylated Antibodies                                                                                                         |
| Application      | WB, ELISA                                                                                                                         |
| Cross Reactivity | Human                                                                                                                             |
| SWISS            | P04626                                                                                                                            |
| Protein Weight   | 138kDa                                                                                                                            |
| Shipping         | Ice bag                                                                                                                           |