

Anti-Phospho-Histone H1.3-T17/Histone H1.4-T17 Monoclonal Antibody For Wb, Elisa - P16402 / P10412

Item Number: CM0001370



Introduction

Rabbit monoclonal antibody specific to phosphorylated T17 of Histone H1.3 and H1.4. Validated for WB and ELISA in human, mouse, and rat samples. Ideal for chromatin research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H1.3-T17/Histone H1.4-T17 Rabbit mAb
Remarks/Alias	H1.3; H1D; H1F3; H1s-2; Phospho-Histone H1.3-T17/Histone H1.4-T17
Species	Human
GeneID	3007/3008
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16402 / P10412
Protein Weight	30 kDa
Shipping	Ice bag