


Anti-Phospho-Histone H2A-S129 Rabbit Pab - P04911

Item Number: CM0000860



Introduction

Rabbit polyclonal antibody targeting phosphorylated Ser129 of Histone H2A, validated for IHC-P and ELISA with human/rat cross-reactivity. Useful for DNA repair, chromatin, and histone code research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H2A-S129 Rabbit pAb
Remarks/Alias	H2A1, SPT11, Phospho-Histone H2A-S129
Species	Human
GeneID (Human)	H2A
GeneID	851811
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S129 of Saccharomyces cerevisiae Histone H2A (NP_010511.3?)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P04911
Protein Weight	14kDa
Shipping	Ice bag