


Anti-Phospho-Histone H2Ax-S139 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16104

Item Number: CM0001445



Introduction

Phospho-Histone H2AX-S139 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat; 15 kDa. For DNA damage and repair pathway analysis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H2AX-S139 Rabbit mAb
Remarks/Alias	H2A.X; H2A/X; H2AFX; Phospho-Histone H2AX-S139
Species	Human
Gene ID (Human)	H2A
Gene ID	3014
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16104
Protein Weight	15 kDa
Shipping	Ice bag