


Anti-Phospho-Histone H3-S10 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16695 / P68431

Item Number: CM0000443



Introduction

Rabbit monoclonal antibody against phosphorylated Histone H3 at Ser10, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Core nucleosome component; phosphorylation marks mitotic chromatin condensation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-S10 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-S10
Species	Human
Gene ID (Human)	8350
Gene ID	8290/8350
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S10 of human Histone H3 (P68431).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag