


Anti-Phospho-Histone H3-S10 Polyclonal Antibody For Wb, IHC-P, IF/ICC, Elisa - Q16695 / P68431

Item Number: CM0001830



Introduction

Phospho-Histone H3-S10 Rabbit Polyclonal Antibody validated for Western blot, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat, and other species. Ideal for studying mitotic signaling and histone code research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-S10 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Phospho-Histone H3-S10
Species	Human
Gene ID (Human)	H3
Gene ID	8290/8350
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S10 of human Histone H3 (NP_003484.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag