


Anti-Phospho-Histone H3-S28 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q16695 / P68431

Item Number: CM0001278



Introduction

Phospho-specific rabbit monoclonal antibody targeting Histone H3 phosphorylated at Serine 28. Validated for Western blot, immunoprecipitation, and ELISA; cross-reacts with human, mouse, and rat. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-S28 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Phospho-Histone H3-S28
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S28 of human Histone H3.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag