


Anti-Phospho-Histone H3-T11 Rabbit Monoclonal Antibody For Wb - Q16695 / P68431

Item Number: CM0001734



Introduction

Phospho-Histone H3-T11 Rabbit Monoclonal Antibody validated for WB, IF, IHC, ChIP, and ELISA. Recognizes human, mouse, and rat Histone H3 phosphorylated at Thr11, a key epigenetic marker. Ideal for chromatin and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3-T11 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A;
Species	Human
GeneID	8290/8350
Immunogen	A synthetic phosphorylated peptide around T11 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag