


Anti-Phospho-Histone H3.3 (Ser31) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P84243

Item Number: CM0000586



Introduction

Rabbit monoclonal antibody detecting phosphorylated Histone H3.3 at Ser31. Suitable for WB, IHC-P, ELISA in human, mouse, rat, and other species. Supports epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3.3-S31 Rabbit mAb
Remarks/Alias	H3-3A; H3.3B; H3F3B; BRYLIB2; Phospho-Histone H3.3-S31
Species	Human
Gene ID (Human)	3021
Gene ID	3020
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S31 of human Histone H3.3 (P84243).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P84243
Protein Weight	15kDa
Shipping	Ice bag