


Anti-Phospho-Histone H3.3-T3 Rabbit Monoclonal Antibody For Wb, Elisa - Q16695 / P68431

Item Number: CM0001484



Introduction

Rabbit monoclonal antibody detecting phosphorylated Thr3 of histone H3.3 (UniProt P84243), validated for WB and ELISA. Cross-reacts with human, mouse, rat, and other species. Ideal for epigenetic research and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Histone H3.3-T3 Rabbit mAb
Remarks/Alias	H3F3; H3-3B; H3.3A; H3F3A; BRYLIB1; Phospho-Histone H3.3-T3
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic phosphorylated peptide around T3 of human Histone H3 (P84243)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag