


Anti-Phospho-Hsl-S563 Rabbit Polyclonal Antibody For Wb, Elisa - Q05469

Item Number: CM0000785



Introduction

Phospho-HSL-S563 rabbit polyclonal antibody for WB and ELISA. Detects phosphorylated HSL (Ser563) in human, mouse, rat samples. Ideal for studying lipolysis, hormone-sensitive lipase activation, and metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSL-S563 Rabbit pAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6; Phospho-HSL-S563
Species	Human
GeneID (human)	3991
GeneID	3991
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S563 of rat HSL (XP_006228456.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag