


Anti-Phospho-Hsl-S660 Rabbit Polyclonal Antibody - Q05469

Item Number: CM0000862



Introduction

Phospho-specific rabbit polyclonal antibody targeting HSL phosphorylated at Ser660. Validated for WB and ELISA; reacts with human and mouse. Detects hormone-sensitive lipase involved in lipid metabolism and energy homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSL-S660 Rabbit pAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6; Phospho-HSL-S660
Species	Human
Gene ID (Human)	3991
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S660 of human LIPE (NP_005348.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag