


Anti-Phospho-Hsl-S853 Rabbit Monoclonal Antibody For Wb - Q05469

Item Number: CM0000652



Introduction

Phospho-specific rabbit monoclonal antibody targeting HSL phosphorylated at serine 853. Validated for Western blot and ELISA, with cross-reactivity to mouse. Essential for studying lipolysis and metabolic signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSL-S853 Rabbit mAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6; Phospho-HSL-S853
Species	Human
GeneID (Human)	3991
GeneID	3991
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S853 of human HSL (Q05469).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag