


Anti-Phospho-Hsl (Ser660) Rabbit Polyclonal Antibody For Wb - Q05469

Item Number: CM0000814



Introduction

Phospho-HSL (Ser660) rabbit polyclonal antibody validated for Western blot. Detects hormone-sensitive lipase phosphorylation in human, mouse, and rat samples. Ideal for lipid metabolism, obesity, and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSL-S660 Rabbit pAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6; Phospho-HSL-S660
Species	Human
GeneID (Human)	3991
GeneID	3991
Immunogen	Synthetic Peptide: A synthesized peptide derived from human HSL (Phospho-Ser660)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag