


Anti-Phospho-Hsp27/Hspb1 (S15) Polyclonal Antibody For Wb, Ip, Elisa - P04792

Item Number: CM0001473



Introduction

Phospho-HSP27/HSPB1 (S15) rabbit polyclonal antibody validated for WB, IP, ELISA. Reacts with human HSP27. Ideal for stress response, chaperone, and actin organization research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S15 Rabbit pAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S15
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S15 of human HSP27/HSPB1 (NP_001531.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag