


Anti-Phospho-Hsp27/Hspb1-S78 Rabbit Monoclonal Antibody For Wb And Elisa - P04792

Item Number: CM0000934



Introduction

Phospho-HSP27/HSPB1-S78 Rabbit Monoclonal Antibody validated for WB and ELISA. Specifically detects human HSPB1 phosphorylated at serine 78. Ideal for chaperone, stress response, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S78 Rabbit mAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S78
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	A synthetic phosphorylated peptide around S78 of human HSPB1 (P04792)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag