

Anti-Phospho-Igf1R (Y1280) Rabbit Polyclonal Antibody For Wb - P08069

Item Number: CM0001551



Introduction

Rabbit polyclonal antibody specific for phosphorylated tyrosine 1280 of human IGF1R. Validated for Western blot, ideal for studying IGF1R signaling in cancer and cell growth research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IGF1R-Y1280 Rabbit pAb
Remarks/Alias	IGFR; CD221; IGFIR; JTK13; Phospho-IGF1R-Y1280
Species	Human
Gene ID	3480
Immunogen	Synthetic peptide corresponding to residues surrounding Y1280 of human IGF1R
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P08069
Protein Weight	155kDa
Shipping	Ice bag