

Anti-Phospho-Ikba (Ser32) Rabbit Monoclonal Antibody For Wb, Elisa - P25963

Item Number: CM0001535



Introduction

Phospho-IkBa (Ser32) rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat IkBa. Detects phosphorylation at Ser32, critical for NF-κB signaling pathway activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IkBa-S32 Rabbit mAb
Remarks/Alias	IKBA; MAD-3; NFKB1; EDAID2; Phospho-IkBa-S32
Species	Human
GeneID (Human)	4792
GeneID	4792
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S32 of human IKB alpha (P25963).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag