

Anti-Phospho-IkB α (Ser36) Rabbit Polyclonal Antibody For Wb And Elisa - P25963

Item Number: CM0001259



Introduction

Phospho-specific rabbit polyclonal antibody targeting human IkB α phosphorylated at Ser36 for Western blot and ELISA applications. Cross-reacts with mouse and rat. Suitable for NF- κ B pathway analysis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IkB α -S36 Rabbit pAb
Remarks/Alias	IKBA; MAD-3; NFKB1; EDAID2; Phospho-IkB α -S36
Species	Human
GeneID (Human)	4792
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S36 of human IkB α (NP_065390.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag