

Anti-Phospho-Insr-Y1150/Igf1R-Y1135 Rabbit Polyclonal Antibody For Wb, Elisa - P06213 / P08069

Item Number: CM0000440



Introduction

Rabbit polyclonal antibody targeting phosphorylated INSR (Y1150) and IGF1R (Y1135), suitable for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for insulin signaling, metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-INSR-Y1150/IGF1R-Y1135 Rabbit pAb
Remarks/Alias	IGF1R/INSR; Phospho-INSR-Y1150/IGF1R-Y1135
Species	Human
GeneID (Human)	3643/3480
GeneID	3643/3480
Immunogen	A synthetic phosphorylated peptide around Y1150 of human IGF1R/INSR (NP_001073285.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06213/P08069
Protein Weight	154kDa/155kDa/156kDa
Shipping	Ice bag