


Anti-Phospho-Irak4 (Thr345/Ser346) Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nwz3

Item Number: CM0001321



Introduction

Specific rabbit polyclonal antibody targeting phosphorylated IRAK4 at Thr345 and Ser346. Validated for Western blotting and ELISA with human samples. Ideal for innate immunity and TLR signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRAK4-Thr345/Ser346 Rabbit pAb
Remarks/Alias	IPD1; IMD67; REN64; IRAK-4; NY-REN-64; Phospho-IRAK4-Thr345/Ser346
Species	Human
Gene ID (Human)	51135
Gene ID	51135
Immunogen	A synthetic phosphorylated peptide around Thr345/Ser346 of human IRAK4.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NWZ3
Protein Weight	52kDa
Shipping	Ice bag