

Anti-Phospho-Irf3-S396 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14653

Item Number: CM0001995



Introduction

Phospho-IRF3-S396 Rabbit Monoclonal Antibody detects phosphorylated IRF3 at Ser396 in human and mouse samples. Validated for WB, IF/ICC, ELISA. Ideal for studying innate immunity, antiviral response, and interferon signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRF3-S396 Rabbit mAb
Remarks/Alias	IIAE7; Phospho-IRF3-S396
Species	Human
Gene ID (Human)	3661
Gene ID	3661
Immunogen	A synthetic phosphorylated peptide around S396 of human IRF3 (NP_001562.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag