

Anti-Phospho-Irf3 (Ser386) Rabbit Monoclonal Antibody For Wb, Elisa - Q14653

Item Number: CM0000573



Introduction

High-specificity rabbit monoclonal antibody against phospho-IRF3 (Ser386) for Western blot and ELISA applications. Detects human IRF3 phosphorylated at Ser386, a key regulator of antiviral innate immune responses.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRF3-S386 Rabbit mAb
Remarks/Alias	IIAE7; Phospho-IRF3-S386
Species	Human
Gene ID (Human)	3661
Gene ID	3661
Immunogen	A synthetic phosphorylated peptide around S386 of human IRF3 (Q14653).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag