

Anti-Phospho-Irs1-S612 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35568

Item Number: CM0000646



Introduction

Phospho-IRS1-S612 Rabbit pAb specifically detects IRS1 phosphorylated at serine 612. Validated for WB, IHC-P, ELISA in human, mouse, rat. Essential tool for insulin signaling and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-IRS1-S612 Rabbit pAb
Remarks/Alias	HIRS-1; Phospho-IRS1-S612
Species	Human
GeneID (Human)	3667
GeneID	3667
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S612 of human IRS1 (NP_005535.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag