

## Anti-Phospho-Irs1-S636 Rabbit Polyclonal Antibody For Wb - P35568

Item Number: CM0000490



### Introduction

Phospho-IRS1-S636 Rabbit Polyclonal Antibody validated for Western blot in human and mouse. Targets IRS1, a central insulin signaling adapter; phosphorylation at serine 636 modulates pathway activity. Suitable for diabetes and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-IRS1-S636 Rabbit pAb                                                                            |
| Remarks/Alias    | HIRS-1; Phospho-IRS1-S636                                                                               |
| Species          | Human                                                                                                   |
| GeneID           | 3667                                                                                                    |
| Immunogen        | A phospho specific peptide corresponding to residues surrounding S636 of human IRS1 (Synthetic peptide) |
| Source           | Rabbit                                                                                                  |
| Category         | Phosphorylated Antibodies                                                                               |
| Application      | WB                                                                                                      |
| Cross Reactivity | Human, Mouse                                                                                            |
| SWISS            | P35568                                                                                                  |
| Protein Weight   | 132kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                 |