

Anti-Phospho-Jak2 (Y1007) Rabbit Polyclonal Antibody For Wb And Elisa - O60674

Item Number: CM0001102



Introduction

Rabbit polyclonal antibody recognizing phosphorylated JAK2 at tyrosine 1007, validated for WB and ELISA in human and mouse. Suitable for studying JAK/STAT signaling, cytokine pathways, and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JAK2-Y1007 Rabbit pAb
Remarks/Alias	JTK10; Phospho-JAK2-Y1007
Species	Human
GeneID (Human)	3717
GeneID	3717
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y1007 of human JAK2 (NP_004963.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag