

Anti-Phospho-Jak2-Y1007 Rabbit Polyclonal Antibody - O60674

Item Number: CM0000991



Introduction

Rabbit polyclonal antibody specific for JAK2 phosphorylated at Tyr1007. Validated for Western blot (WB) and ELISA, reacts with human and mouse JAK2. Ideal for studying JAK-STAT signaling, cytokine responses, and chromatin regulation.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-JAK2-Y1007 Rabbit pAb
Remarks/Alias	JTK10; Phospho-JAK2-Y1007
Species	Human
GeneID (Human)	3717
GeneID	3717
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around Y1007 of human JAK2 (NP_004963.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag