

Anti-Phospho-Jak2 (Y1007/1008) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60674

Item Number: CM0001210



Introduction

Rabbit polyclonal antibody specific for phosphorylated JAK2 (Y1007/1008). Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Suited for cytokine signaling and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JAK2-Y1007/1008 Rabbit pAb
Remarks/Alias	JTK10; Phospho-JAK2-Y1007/1008
Species	Human
GeneID (human)	3717
GeneID	3717
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around Y1007 & Y1008 of human JAK2 (NP_004963.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60674
Protein Weight	131kDa
Shipping	Ice bag