

Anti-Phospho-Jnk1-T183/Y185 + Jnk2-T183/Y185 + Jnk3-T221/Y223 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P45983/P45984/P53779

Item Number: CM0001391



Introduction

Phospho-specific rabbit monoclonal antibody recognizing activated JNK1/2/3 (Thr183/Tyr185 and Thr221/Tyr223). Suitable for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for stress and apoptosis signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-JNK1-T183/Y185 + JNK2-T183/Y185 + JNK3-T221/Y223 Rabbit mAb
Remarks/Alias	JNK1/2/3; SAPK; Phospho-JNK1-T183/Y185 + JNK2-T183/Y185 + JNK3-T221/Y223
Species	Human
GenelD (Human)	5599/5601/5602
GenelD	5599/5601/5602
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45983/P45984/P53779
Protein Weight	48kDa/52kDa
Shipping	Ice bag