

## Anti-Phospho-Jnk1/2/3 (T183/T183/T221) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P45983 / P45984 / P53779

Item Number: CM0001380



### Introduction

Phospho-specific rabbit monoclonal antibody recognizing JNK1/2/3 phosphorylated at T183/T183/T221. Suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying JNK-mediated stress signaling, apoptosis, and inflammation.

[Learn More](#)

| Parameter        | Value                                                                  |
|------------------|------------------------------------------------------------------------|
| Product Name     | Phospho-JNK1/2/3-T183/T183/T221 Rabbit mAb                             |
| Remarks/Alias    | JNK1/JNK2/JNK3; Phospho-JNK1/2/3-T183/T183/T221                        |
| Species          | Human                                                                  |
| Gene ID (Human)  | 5599/5601/5602                                                         |
| Immunogen        | A synthetic phosphorylated peptide around T183 of human JNK1 (P45983). |
| Source           | Rabbit                                                                 |
| Category         | Phosphorylated Antibodies                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                               |
| Cross Reactivity | Human, Mouse, Rat                                                      |
| SWISS            | P45983/P45984/P53779                                                   |
| Protein Weight   | 35kDa/44kDa/48kDa/27kDa/52kDa                                          |
| Shipping         | Ice bag                                                                |